


## Individual Zircon Abutment with Titanbase.



**Titanbase**






Implant name	Implant type	Diameter	Titanbasis Code
Nobel Biocare	Replace select NP	3.5 mm	E 800
	Replace select RP	4.3 mm	E 810
	Replace select WP	5.0 mm	E 820
	Replace select WP	6.0 mm	E 830
Biomet 3i	Osseotite	3.4 mm	I 800
	Osseotite	4.1 mm	IK 810
	Osseotite	5.0 mm	I 820
Nobel Biocare Brånemark	Brånemark NP	3.4 mm	K 800
	Brånemark RP	4,1 mm	IK 810
	Brånemark WP	5.1 mm	K 820
Straumann Implantat	SynOcta NN	3.5 mm	N 800
	SynOcta RN	4.8 mm	N 810
	SynOcta WN	6.5 mm	N 820
Astra	Osseospeed Aqua 3.5/4.0	3.5/4.0 mm	S 800
	Osseospeed Lila 4.5/5.0	4.5/5.0 mm	S 820



Mouth Scanner

# iTero Mouth Scanner





The coloured areas are distance to antagonist less as 1mm



Milled model



First case from: Itero scanner /  
Bite Europe / Lotus dental / Hexa  
ceram





Sent the data to USA /German  
in future to Lotus dental for  
milling a model.

Server to download  
or upload data from  
the dentist and  
Cadent.

Cadent itero data is  
available for 3Shape  
Dental System  
Software





Zirconium dioxide:  $ZrO_2$  Y-TZP-A  $ZrO_2 + HfO_2 + Y_2O_3$

Caping /bridges /Veneer /Inlay / Maryland /Abutment  
up to 14 units bridge

In 4 different color / **\*NEW\*** Zirconium dioxide with full colour



Polyurethane / Wax crowns and bridges for casting  
up to 14 units bridge



Temporary crowns and bridges / PMMA

up to 14 units bridge

In A2 / A3 /B1 / B2 color



**\*NEW\*** Individual Zircon Abutment with Titanbase.

Zirconium dioxide:  $ZrO_2$  Y-TZP-A  $ZrO_2 + HfO_2 + Y_2O_3$

Titanbase Medenika Germany