

## Lotus Individual Zircon Abutment with Titanbase 2012

Implant type	Product Code	Description	Diameter
Nobel Biocare Replace	Abutment:E 810/Screw:E 61	Medentika GmbH	4.3 mm.
Nobel Biocare Replace	Abutment:E 800/Screw:E 60	Medentika GmbH	3.5 mm.
Nobel Biocare Replace	Abutment:E 820/Screw:E 61	Medentika GmbH	5.0 mm.
Nobel Biocare Replace	Abutment:E 830/Screw:E 61	Medentika GmbH	6.0 mm.
Biomet 3i	Abutment:H 800/Screw:H 60	Medentika GmbH	3.4 mm.
Biomet 3i	Abutment:H 810/Screw:H 60	Medentika GmbH	4.1 mm.
Biomet 3i	Abutment:H 820/Screw:H 60	Medentika GmbH	5.0 mm.
Biomet 3i External Hex	Abutment:I 800/Screw:I 61	Medentika GmbH	3.4 mm.
Biomet 3i External Hex	Abutment:I 810/Screw:I 61	Medentika GmbH	4.1 mm.
Biomet 3i External Hex	Abutment:I 820/Screw:I 61	Medentika GmbH	5.0 mm.
Nobel Biocare Brancmark	Abutment:K 800/Screw:K 60	Medentika GmbH	3.5 mm.
Nobel Biocare Brancmark	Abutment:K 810/Screw:K 61	Medentika GmbH	4.1 mm.
Nobel Biocare Brancmark	Abutment:K 820/Screw:K 62	Medentika GmbH	5.1 mm.
Straumann Bone Level	Abutment:L 810/Screw:L 60	Medentika GmbH	
Straumann Bone Level	Abutment:L 811/Screw:L 60	Medentika GmbH	
Straumann Bone Level	Abutment:L 800/Screw:L 60	Medentika GmbH	
Straumann Bone Level	Abutment:L 801/Screw:L 60	Medentika GmbH	
Stauamm synocta	Abutment:N 800/Screw:N 60	Medentika GmbH	3.5 mm.
Stauamm synocta	Abutment:N 810/Screw:N 62	Medentika GmbH	4.8 mm.
Stauamm synocta	Abutment:N 820/Screw:N 62	Medentika GmbH	6.5 mm.
Zimmer Tapered Screw - Vent	Abutment:R 810/Screw:R 60	Medentika GmbH	4.5 mm.
Zimmer Tapered Screw - Vent	Abutment:R 800/Screw:R 60	Medentika GmbH	3.5 mm.
Zimmer Tapered Screw - Vent	Abutment:R 820/Screw:R 60	Medentika GmbH	5.7 mm.
Astra Tech Osseo Speed	Abutment:S 820/Screw:S 61	Medentika GmbH	4.5-5.0 mm.
Astra Tech Osseo Speed	Abutment:S 800/Screw:S 60	Medentika GmbH	3.5-4.0 mm.
Dentsply Friadent Frialit / Xive	Abutment:T 800/Screw:T 60	Medentika GmbH	-
Dentsply Friadent Frialit / Xive	Abutment:T 805/Screw:T 60	Medentika GmbH	3.8 mm.
Dentsply Friadent Frialit / Xive	Abutment:T 810/Screw:T 60	Medentika GmbH	-
Dentsply Friadent Frialit / Xive	Abutment:T 820/Screw:T 60	Medentika GmbH	5.5 mm.
Nobel Active	Abutment: F 800/Screw:F 60	Medentika GmbH	Torque : 25 ncm
Nobel Active	Abutment: F810/Screw:F 61	Medentika GmbH	Torque : 35 ncm
Camlog	Abutment: Titanium basis with Screw : 3.3 mm.	Camlog	3.3 mm.
Camlog	Abutment: Titanium basis with Screw : 3.8 mm.	Camlog	3.8 mm.
Camlog	Abutment: Titanium basis with Screw : 4.3 mm.	Camlog	4.3 mm.
Camlog	Abutment: Titanium basis with Screw : 5.0 mm.	Camlog	5.0 mm.
Camlog	Abutment: Titanium basis with Screw : 6.0 mm.	Camlog	6.0 mm.
THOMMEN	Abutment screw included 4.0	THOMMEN	4.0 mm.
THOMMEN	Abutment screw included 4.5	THOMMEN	4.5 mm.
THOMMEN	Abutment screw included 5.0	THOMMEN	5.0 mm.
Biocomp		Biocomp	
Ankylos	C standard	Ankylos	
Ankylos	C high	Ankylos	
Ankylos	X Standard	Ankylos	
Ankylos	X high	Ankylos	
MegaGen AnyRidge Internal	AnyRidge	MegaGen	4.0-5.5mm
MegaGen EZ Plus Internal	Regular Type	MegaGen	
MegaGen EZ Plus Internal	Small Type	MegaGen	
MegaGen Rescue Internal	Rescue	MegaGen	6.0mm

**NEW !**